

For Research Use Only

# NINJ2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50693-1

## Capture Antibody Information

Catalog Number:  
66202-2-PBS

Host:  
Mouse

Isotype:  
IgG2b

Purification Method:  
Protein A Magarose purification

Clone ID:  
2E4B1

Reactivity:  
human

Immunogen Catalog Number:  
Ag5343

Conjugate:  
Unconjugated

Full name:  
ninjurin 2

Gene ID:  
4815

## Detection Antibody Information

Catalog Number:  
66202-3-PBS

Host:  
Mouse

Isotype:  
IgG2b

Purification Method:  
Protein A purification

Clone ID:  
1D1B10

Reactivity:  
human

GenBank:  
BC057766

Immunogen Catalog Number:  
Ag5343

Conjugate:  
Unconjugated

Full name:  
ninjurin 2

Gene ID:  
4815

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50693-1 targets NINJ2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NINJ2 Monoclonal antibody, PBS Only (Capture) 66202-2-PBS (2E4B1). 100 µg. Concentration 1 mg/mL.

Detection antibody: NINJ2 Monoclonal antibody, PBS Only (Detector) 66202-3-PBS (1D1B10). 100 µg. Concentration 1 mg/mL.

Alternative NINJ2 matched antibody pairs: MP50693-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50693-1, NINJ2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66202-2-PBS. Detection antibody: 66202-3-PBS. Standard: Ag5343. Range: 0.781-100 ng/mL.