

NINJ2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50693-2

Capture Antibody Information

Catalog Number:

66202-4-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G Magarose purification

Clone ID:

5E8C2

Reactivity:

human

Immunogen Catalog Number:

Ag5343

Conjugate:

Unconjugated

Full name:

ninjurin 2

Gene ID:

4815

Detection Antibody Information

Catalog Number:

66202-3-PBS

Host:

Mouse

Isotype:

IgG2b

Purification Method:

Protein A purification

Clone ID:

1D1B10

Reactivity:

human

GenBank:

BC057766

Immunogen Catalog Number:

Ag5343

Conjugate:

Unconjugated

Full name:

ninjurin 2

Gene ID:

4815

Applications

Tested Applications:

Cytometric bead array

Range:

0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50693-2 targets NINJ2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NINJ2 Monoclonal antibody, PBS Only (Capture) 66202-4-PBS (5E8C2). 100 µg. Concentration 1 mg/ml.

Detection antibody: NINJ2 Monoclonal antibody, PBS Only (Detector) 66202-3-PBS (1D1B10). 100 µg. Concentration 1 mg/ml.

Alternative NINJ2 matched antibody pairs: MP50693-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50693-2, NINJ2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66202-4-PBS. Detection antibody: 66202-3-PBS. Standard: Ag5343. Range: 0.391-100 ng/mL.