

For Research Use Only

CD44 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50724-3

Capture Antibody Information

Catalog Number:
60224-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
1G7G4

Reactivity:
human

Immunogen Catalog Number:
Ag7633

Conjugate:
Unconjugated

Full name:
CD44 molecule (Indian blood group)

Gene ID:
960

Detection Antibody Information

Catalog Number:
60224-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
8A3F6

Reactivity:
human

GenBank:
BC004372

Immunogen Catalog Number:
Ag7633

Conjugate:
Unconjugated

Full name:
CD44 molecule (Indian blood group)

Gene ID:
960

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50724-3 targets CD44 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CD44 Monoclonal antibody, PBS Only (Capture) 60224-2-PBS (1G7G4). 100 µg. Concentration 1 mg/ml.

Detection antibody: CD44 Monoclonal antibody, PBS Only (Capture/Detector) 60224-4-PBS (8A3F6). 100 µg. Concentration 1 mg/ml.

Alternative CD44 matched antibody pairs: MP00033-1, MP00033-2, MP01236-1, MP01236-2, MP01236-3, MP50724-1, MP50724-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50724-3, CD44 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60224-2-PBS. Detection antibody: 60224-4-PBS. Standard: Ag7633. Range: 0.098-100 ng/mL.