

For Research Use Only

BCL7B Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50729-4

Capture Antibody Information

Catalog Number:
60517-5-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1D1H8

Reactivity:
human

Immunogen Catalog Number:
Ag22892

Conjugate:
Unconjugated

Full name:
B-cell CLL/lymphoma 7B

Gene ID:
9275

Detection Antibody Information

Catalog Number:
60517-2-PBS

Host:
Mouse

Isotype:
IgG2b

Purification Method:
Protein A Magarose purification

Clone ID:
3B11F11

Reactivity:
human

GenBank:
BC000956

Immunogen Catalog Number:
Ag22892

Conjugate:
Unconjugated

Full name:
B-cell CLL/lymphoma 7B

Gene ID:
9275

Applications

Tested Applications:
Cytometric bead array

Range:
0.195-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50729-4 targets BCL7B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BCL7B Monoclonal antibody, PBS Only (Capture) 60517-5-PBS (1D1H8). 100 µg. Concentration 1 mg/ml.

Detection antibody: BCL7B Monoclonal antibody, PBS Only (Detector) 60517-2-PBS (3B11F11). 100 µg. Concentration 1 mg/ml.

Alternative BCL7B matched antibody pairs: MP50729-1, MP50729-2, MP50729-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50729-4, BCL7B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60517-5-PBS. Detection antibody: 60517-2-PBS. Standard: Ag22892. Range: 0.195-100 ng/mL.