

For Research Use Only

GRIM19 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50735-1

Capture Antibody Information

Catalog Number:
67414-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
3C11D4

Reactivity:
human

Immunogen Catalog Number:
Ag20855

Conjugate:
Unconjugated

Full name:
NADH dehydrogenase (ubiquinone) 1
alpha subcomplex, 13

Gene ID:
51079

Detection Antibody Information

Catalog Number:
67414-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
3A10C5

Reactivity:
human

GenBank:
BC009189

Immunogen Catalog Number:
Ag20855

Conjugate:
Unconjugated

Full name:
NADH dehydrogenase (ubiquinone) 1
alpha subcomplex, 13

Gene ID:
51079

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP50735-1 targets GRIM19 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GRIM19 Monoclonal antibody, PBS Only (Capture) 67414-2-PBS (3C11D4). 100 µg. Concentration 1 mg/ml.

Detection antibody: GRIM19 Monoclonal antibody, PBS Only (Capture/Detector) 67414-1-PBS (3A10C5). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50735-1, GRIM19 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67414-2-PBS. Detection antibody: 67414-1-PBS. Standard: Ag20855. Range: 0.098-100 ng/mL.