

For Research Use Only

PDGFRB Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50756-2

Capture Antibody Information

Catalog Number:
60535-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2E6H12

Reactivity:
human

Conjugate:
Unconjugated

Full name:
platelet-derived growth factor receptor, beta polypeptide

Gene ID:
5159

Detection Antibody Information

Catalog Number:
60535-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
3G3G10

Reactivity:
human

GenBank:
BC032224

Conjugate:
Unconjugated

Full name:
platelet-derived growth factor receptor, beta polypeptide

Gene ID:
5159

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50756-2 targets PDGFRB in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PDGFRB Monoclonal antibody, PBS Only (Capture) 60535-3-PBS (2E6H12). 100 µg. Concentration 1 mg/ml.

Detection antibody: PDGFRB Monoclonal antibody, PBS Only (Detector) 60535-4-PBS (3G3G10). 100 µg. Concentration 1 mg/ml.

Alternative PDGFRB matched antibody pairs: MP00006-2, MP00006-3, MP50756-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50756-2, PDGFRB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60535-3-PBS. Detection antibody: 60535-4-PBS. Standard:Eg0147. Range: 0.098-100 ng/mL.