

# VAPA Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50773-1

## Capture Antibody Information

**Catalog Number:**

67845-2-PBS

**Host:**

Mouse

**Isotype:**

IgG1

**Purification Method:**

Protein G Magarose purification

**Clone ID:**

1D2D8

**Reactivity:**

human

**Immunogen Catalog Number:**

Ag7819

**Conjugate:**

Unconjugated

**Full name:**

VAMP (vesicle-associated membrane protein)-associated protein A, 33kDa

**Gene ID:**

9218

## Detection Antibody Information

**Catalog Number:**

67845-3-PBS

**Host:**

Mouse

**Isotype:**

IgG1

**Purification Method:**

Protein G Magarose purification

**Clone ID:**

1H1E6

**Reactivity:**

human

**GenBank:**

BC002992

**Immunogen Catalog Number:**

Ag7819

**Conjugate:**

Unconjugated

**Full name:**

VAMP (vesicle-associated membrane protein)-associated protein A, 33kDa

**Gene ID:**

9218

## Applications

**Tested Applications:**

Cytometric bead array

**Range:**

1.563-200 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50773-1 targets VAPA in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: VAPA Monoclonal antibody, PBS Only (Capture) 67845-2-PBS (1D2D8). 100 µg. Concentration 1 mg/ml.

Detection antibody: VAPA Monoclonal antibody, PBS Only (Detector) 67845-3-PBS (1H1E6). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**

Store at -80°C.

**Storage buffer:**

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50773-1, VAPA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67845-2-PBS. Detection antibody: 67845-3-PBS. Standard: Ag7819. Range: 1.563-200 ng/mL.