

For Research Use Only

HAX1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50775-2

Capture Antibody Information

Catalog Number:
67883-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
2A12F11

Reactivity:
human, mouse, rat

Immunogen Catalog Number:
Ag27244

Conjugate:
Unconjugated

Full name:
HCLS1 associated protein X-1

Gene ID:
10456

Detection Antibody Information

Catalog Number:
67883-5-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2A6F7

Reactivity:
human

GenBank:
BC014314

Immunogen Catalog Number:
Ag27244

Conjugate:
Unconjugated

Full name:
HCLS1 associated protein X-1

Gene ID:
10456

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50775-2 targets HAX1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HAX1 Monoclonal antibody, PBS Only (Capture) 67883-2-PBS (2A12F11). 100 µg. Concentration 1 mg/mL.

Detection antibody: HAX1 Monoclonal antibody, PBS Only (Detector) 67883-5-PBS (2A6F7). 100 µg. Concentration 1 mg/mL.

Alternative HAX1 matched antibody pairs: MP50775-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50775-2, HAX1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67883-2-PBS. Detection antibody: 67883-5-PBS. Standard: Ag27244. Range: 0.391-50 ng/mL