

For Research Use Only

# HAX1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50775-3

## Capture Antibody Information

Catalog Number:  
67883-6-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G purification

Clone ID:  
1F8H3

Reactivity:  
human

Immunogen Catalog Number:  
Ag27244

Conjugate:  
Unconjugated

Full name:  
HCLS1 associated protein X-1

Gene ID:  
10456

## Detection Antibody Information

Catalog Number:  
67883-7-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
2E11H5

Reactivity:  
human

GenBank:  
BC014314

Immunogen Catalog Number:  
Ag27244

Conjugate:  
Unconjugated

Full name:  
HCLS1 associated protein X-1

Gene ID:  
10456

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50775-3 targets HAX1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HAX1 Monoclonal antibody, PBS Only (Capture) 67883-6-PBS (1F8H3). 100 µg. Concentration 1 mg/ml.

Detection antibody: HAX1 Monoclonal antibody, PBS Only (Detector) 67883-7-PBS (2E11H5). 100 µg. Concentration 1 mg/ml.

Alternative HAX1 matched antibody pairs: MP50775-1, MP50775-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50775-3, HAX1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67883-6-PBS. Detection antibody: 67883-7-PBS. Standard: Ag27244. Range: 0.781-100 ng/mL.