

For Research Use Only

TET2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50789-4

Capture Antibody Information

Catalog Number:
60560-7-PBS

Host:
Mouse

Isotype:
IgG2b

Purification Method:
Protein A Magarose purification

Clone ID:
1H8F3

Reactivity:
human

Immunogen Catalog Number:
Ag15773

Conjugate:
Unconjugated

Full name:
tet oncogene family member 2

Gene ID:
54790

Detection Antibody Information

Catalog Number:
60560-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2E5A8

Reactivity:
human

GenBank:
BC150180

Immunogen Catalog Number:
Ag15773

Conjugate:
Unconjugated

Full name:
tet oncogene family member 2

Gene ID:
54790

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50789-4 targets TET2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TET2 Monoclonal antibody, PBS Only (Capture) 60560-7-PBS (1H8F3). 100 µg. Concentration 1 mg/ml.

Detection antibody: TET2 Monoclonal antibody, PBS Only (Detector) 60560-2-PBS (2E5A8). 100 µg. Concentration 1 mg/ml.

Alternative TET2 matched antibody pairs: MP50789-1, MP50789-2, MP50789-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50789-4, TET2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60560-7-PBS. Detection antibody: 60560-2-PBS. Standard: Ag15773. Range: 0.391-100 ng/mL.