

For Research Use Only

CREB3L1,OASIS Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50829-3

Capture Antibody Information

Catalog Number:
67617-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1G6A2
Reactivity:
human
Immunogen Catalog Number:
Ag28042

Conjugate:
Unconjugated
Full name:
cAMP responsive element binding protein 3-like 1
Gene ID:
90993

Detection Antibody Information

Catalog Number:
67617-5-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1B4A9
Reactivity:
human
GenBank:
BC015781
Immunogen Catalog Number:
Ag28042

Conjugate:
Unconjugated
Full name:
cAMP responsive element binding protein 3-like 1
Gene ID:
90993

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50829-3 targets CREB3L1,OASIS in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CREB3L1,OASIS Monoclonal antibody, PBS Only (Capture) 67617-2-PBS (1G6A2). 100 µg. Concentration 1 mg/ml.

Detection antibody: CREB3L1,OASIS Monoclonal antibody, PBS Only (Detector) 67617-5-PBS (1B4A9). 100 µg. Concentration 1 mg/ml.

Alternative CREB3L1,OASIS matched antibody pairs: MP50829-1, MP50829-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50829-3, CREB3L1, OASIS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67617-2-PBS. Detection antibody: 67617-5-PBS. Standard: Ag28042. Range: 1.563-100 ng/mL.