

For Research Use Only

HSD17B2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50836-3

Capture Antibody Information

Catalog Number:
60587-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2C1A2
Reactivity:
human
Immunogen Catalog Number:
Ag1411

Conjugate:
Unconjugated
Full name:
hydroxysteroid (17-beta) dehydrogenase 2
Gene ID:
3294

Detection Antibody Information

Catalog Number:
60587-4-PBS
Host:
Mouse
Isotype:
IgG2b
Purification Method:
Protein A purification

Clone ID:
2B10B4
Reactivity:
human
GenBank:
BC009581
Immunogen Catalog Number:
Ag1411

Conjugate:
Unconjugated
Full name:
hydroxysteroid (17-beta) dehydrogenase 2
Gene ID:
3294

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50836-3 targets HSD17B2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: HSD17B2 Monoclonal antibody, PBS Only (Capture) 60587-1-PBS (2C1A2). 100 µg. Concentration 1 mg/ml.
Detection antibody: HSD17B2 Monoclonal antibody, PBS Only (Detector) 60587-4-PBS (2B10B4). 100 µg. Concentration 1 mg/ml.
Alternative HSD17B2 matched antibody pairs: MP50836-1, MP50836-2
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50836-3, HSD17B2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60587-1-PBS. Detection antibody: 60587-4-PBS. Standard:Ag1411. Range: 0.781-100 ng/mL.