

For Research Use Only

NEDD4 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50851-2

Capture Antibody Information

Catalog Number:
67865-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
3D3A3

Reactivity:
human

Immunogen Catalog Number:
Ag17124

Conjugate:
Unconjugated

Full name:
neural precursor cell expressed,
developmentally down-regulated 4

Gene ID:
4734

Detection Antibody Information

Catalog Number:
67865-5-PBS

Host:
Mouse

Isotype:
IgG3

Purification Method:
Protein A Magarose purification

Clone ID:
1B10H5

Reactivity:
human

GenBank:
BC152452

Immunogen Catalog Number:
Ag17124

Conjugate:
Unconjugated

Full name:
neural precursor cell expressed,
developmentally down-regulated 4

Gene ID:
4734

Applications

Tested Applications:
Cytometric bead array

Range:
6.25-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50851-2 targets NEDD4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NEDD4 Monoclonal antibody, PBS Only (Capture) 67865-4-PBS (3D3A3). 100 µg. Concentration 1 mg/ml.

Detection antibody: NEDD4 Monoclonal antibody, PBS Only (Detector) 67865-5-PBS (1B10H5). 100 µg. Concentration 1 mg/ml.

Alternative NEDD4 matched antibody pairs: MP00234-1, MP00234-2, MP50851-1, MP50851-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50851-2, NEDD4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67865-4-PBS. Detection antibody: 67865-5-PBS. Standard: Ag17124. Range: 6.25-100 ng/mL