

For Research Use Only

# RPL28 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50863-3

## Capture Antibody Information

Catalog Number:

60607-1-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G Magarose purification

Clone ID:

1F5E1

Reactivity:

human

Immunogen Catalog Number:

Ag9229

Conjugate:

Unconjugated

Full name:

ribosomal protein L28

Gene ID:

6158

## Detection Antibody Information

Catalog Number:

60607-3-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G purification

Clone ID:

1E4B2

Reactivity:

human

GenBank:

BC010173

Immunogen Catalog Number:

Ag9229

Conjugate:

Unconjugated

Full name:

ribosomal protein L28

Gene ID:

6158

## Applications

Tested Applications:

Cytometric bead array

Range:

0.781-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50863-3 targets RPL28 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RPL28 Monoclonal antibody, PBS Only (Capture) 60607-1-PBS (1F5E1). 100 µg. Concentration 1 mg/ml.

Detection antibody: RPL28 Monoclonal antibody, PBS Only (Detector) 60607-3-PBS (1E4B2). 100 µg. Concentration 1 mg/ml.

Alternative RPL28 matched antibody pairs: MP50863-1, MP50863-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50863-3, RPL28 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60607-1-PBS. Detection antibody: 60607-3-PBS. Standard: Ag9229. Range: 0.781-50 ng/mL