

For Research Use Only

# SILV Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50869-3

## Capture Antibody Information

Catalog Number:  
60612-4-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
2A6H6

Reactivity:  
human

Immunogen Catalog Number:  
Ag6726

Conjugate:  
Unconjugated

Full name:  
silver homolog (mouse)

Gene ID:  
6490

## Detection Antibody Information

Catalog Number:  
60612-5-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
1C10F8

Reactivity:  
human

GenBank:  
BC001414

Immunogen Catalog Number:  
Ag6726

Conjugate:  
Unconjugated

Full name:  
silver homolog (mouse)

Gene ID:  
6490

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.098-6.25 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50869-3 targets SILV in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SILV Monoclonal antibody, PBS Only (Capture) 60612-4-PBS (2A6H6). 100 µg. Concentration 1 mg/ml.

Detection antibody: SILV Monoclonal antibody, PBS Only (Detector) 60612-5-PBS (1C10F8). 100 µg. Concentration 1 mg/ml.

Alternative SILV matched antibody pairs: MP50869-1, MP50869-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50869-3, SILV Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60612-4-PBS. Detection antibody: 60612-5-PBS. Standard:Ag6726. Range: 0.098-6.25 ng/mL