

For Research Use Only

TEX264 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50870-2

Capture Antibody Information

Catalog Number:

68507-2-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G Magarose purification

Clone ID:

1B5B10

Reactivity:

human

Immunogen Catalog Number:

Ag33710

Conjugate:

Unconjugated

Full name:

testis expressed 264

Gene ID:

51368

Detection Antibody Information

Catalog Number:

68507-3-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G purification

Clone ID:

1D8H8

Reactivity:

human

GenBank:

BC008742

Immunogen Catalog Number:

Ag33710

Conjugate:

Unconjugated

Full name:

testis expressed 264

Gene ID:

51368

Applications

Tested Applications:

Cytometric bead array

Range:

0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50870-2 targets TEX264 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TEX264 Monoclonal antibody, PBS Only (Capture/Detector) 68507-2-PBS (1B5B10). 100 µg. Concentration 1 mg/ml.

Detection antibody: TEX264 Monoclonal antibody, PBS Only (Detector) 68507-3-PBS (1D8H8). 100 µg. Concentration 1 mg/ml.

Alternative TEX264 matched antibody pairs: MP50870-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50870-2, TEX264 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68507-2-PBS. Detection antibody: 68507-3-PBS. Standard: Ag33710. Range: 0.098-100 ng/mL.