

For Research Use Only

KIAA0776 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50874-1

Capture Antibody Information

Catalog Number:

60614-1-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G Magarose purification

Clone ID:

3E9G8

Reactivity:

human

Immunogen Catalog Number:

Ag23404

Conjugate:

Unconjugated

Full name:

KIAA0776

Gene ID:

23376

Detection Antibody Information

Catalog Number:

60614-2-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G purification

Clone ID:

3G11F6

Reactivity:

human

GenBank:

BC036379

Immunogen Catalog Number:

Ag23404

Conjugate:

Unconjugated

Full name:

KIAA0776

Gene ID:

23376

Applications

Tested Applications:

Cytometric bead array

Range:

0.195-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50874-1 targets KIAA0776 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: KIAA0776 Monoclonal antibody, PBS Only (Capture) 60614-1-PBS (3E9G8). 100 µg. Concentration 1 mg/ml.

Detection antibody: KIAA0776 Monoclonal antibody, PBS Only (Detector) 60614-2-PBS (3G11F6). 100 µg. Concentration 1 mg/ml.

Alternative KIAA0776 matched antibody pairs: MP00135-1, MP00135-2, MP00135-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50874-1, KIAA0776 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60614-1-PBS. Detection antibody: 60614-2-PBS. Standard: Ag23404. Range: 0.195-200 ng/mL