

For Research Use Only

# SEC11A Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50880-3

## Capture Antibody Information

Catalog Number:  
67379-5-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
1F9D3

Reactivity:  
human

Immunogen Catalog Number:  
Ag6610

Conjugate:  
Unconjugated

Full name:  
SEC11 homolog A (*S. cerevisiae*)

Gene ID:  
23478

## Detection Antibody Information

Catalog Number:  
67379-6-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
2E1E2

Reactivity:  
human

GenBank:  
BC000359

Immunogen Catalog Number:  
Ag6610

Conjugate:  
Unconjugated

Full name:  
SEC11 homolog A (*S. cerevisiae*)

Gene ID:  
23478

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.098-12.5 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50880-3 targets SEC11A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SEC11A Monoclonal antibody, PBS Only (Capture) 67379-5-PBS (1F9D3). 100 µg. Concentration 1 mg/ml.

Detection antibody: SEC11A Monoclonal antibody, PBS Only (Detector) 67379-6-PBS (2E1E2). 100 µg. Concentration 1 mg/ml.

Alternative SEC11A matched antibody pairs: MP50880-1, MP50880-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50880-3, SEC11A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67379-5-PBS. Detection antibody: 67379-6-PBS. Standard:Ag6610. Range: 0.098-12.5 ng/mL

Cytometric bead array standard curve of MP50880-3, SEC11A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67379-5-PBS. Detection antibody: 67379-6-PBS. Standard:null. Range: null.