

For Research Use Only

GRAMD1B Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50885-2

Capture Antibody Information

Catalog Number:
67551-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1D1A3

Reactivity:
human

Immunogen Catalog Number:
Ag20670

Conjugate:
Unconjugated

Full name:
GRAM domain containing 1B

Gene ID:
57476

Detection Antibody Information

Catalog Number:
67551-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1C2F2

Reactivity:
human

GenBank:
BC107480

Immunogen Catalog Number:
Ag20670

Conjugate:
Unconjugated

Full name:
GRAM domain containing 1B

Gene ID:
57476

Applications

Tested Applications:
Cytometric bead array

Range:
0.049-3.125 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50885-2 targets GRAMD1B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GRAMD1B Monoclonal antibody, PBS Only (Capture) 67551-4-PBS (1D1A3). 100 µg. Concentration 1 mg/ml.

Detection antibody: GRAMD1B Monoclonal antibody, PBS Only (Detector) 67551-3-PBS (1C2F2). 100 µg. Concentration 1 mg/ml.

Alternative GRAMD1B matched antibody pairs: MP50885-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50885-2, GRAMD1B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67551-4-PBS. Detection antibody: 67551-3-PBS. Standard: Ag20670. Range: 0.049-3.125 ng/mL.