

For Research Use Only

GRAMD1B Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50885-4

Capture Antibody Information

Catalog Number:
67551-7-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
1H9H2

Reactivity:
human

Immunogen Catalog Number:
Ag20670

Conjugate:
Unconjugated

Full name:
GRAM domain containing 1B

Gene ID:
57476

Detection Antibody Information

Catalog Number:
67551-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1C2F2

Reactivity:
human

GenBank:
BC107480

Immunogen Catalog Number:
Ag20670

Conjugate:
Unconjugated

Full name:
GRAM domain containing 1B

Gene ID:
57476

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50885-4 targets GRAMD1B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GRAMD1B Monoclonal antibody, PBS Only (Capture) 67551-7-PBS (1H9H2). 100 µg. Concentration 1 mg/ml.

Detection antibody: GRAMD1B Monoclonal antibody, PBS Only (Detector) 67551-3-PBS (1C2F2). 100 µg. Concentration 1 mg/ml.

Alternative GRAMD1B matched antibody pairs: MP50885-1, MP50885-2, MP50885-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50885-4, GRAMD1B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67551-7-PBS. Detection antibody: 67551-3-PBS. Standard: Ag20670. Range: 1.563-50 ng/mL