

For Research Use Only

# TRIM36 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50936-3

## Capture Antibody Information

Catalog Number:  
60653-4-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G purification

Clone ID:  
3C7B8

Reactivity:  
human

Immunogen Catalog Number:  
Ag23227

Conjugate:  
Unconjugated

Full name:  
tripartite motif-containing 36

Gene ID:  
55521

## Detection Antibody Information

Catalog Number:  
60653-5-PBS

Host:  
Mouse

Isotype:  
IgG3

Purification Method:  
Protein A purification

Clone ID:  
2F8B6

Reactivity:  
human

GenBank:  
BC130334

Immunogen Catalog Number:  
Ag23227

Conjugate:  
Unconjugated

Full name:  
tripartite motif-containing 36

Gene ID:  
55521

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.781-25 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50936-3 targets TRIM36 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TRIM36 Monoclonal antibody, PBS Only (Capture) 60653-4-PBS (3C7B8). 100 µg. Concentration 1 mg/mL.

Detection antibody: TRIM36 Monoclonal antibody, PBS Only (Detector) 60653-5-PBS (2F8B6). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50936-3, TRIM36 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60653-4-PBS. Detection antibody: 60653-5-PBS. Standard: Ag23227. Range: 0.781-25 ng/mL