

For Research Use Only

DGAT1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50948-3

Capture Antibody Information

Catalog Number:
60664-5-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1C9D3
Reactivity:
human
Immunogen Catalog Number:
Ag32708

Conjugate:
Unconjugated
Full name:
diacylglycerol O-acyltransferase homolog 1 (mouse)
Gene ID:
8694

Detection Antibody Information

Catalog Number:
60664-6-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2G5H12
Reactivity:
human
GenBank:
BC015762
Immunogen Catalog Number:
Ag32708

Conjugate:
Unconjugated
Full name:
diacylglycerol O-acyltransferase homolog 1 (mouse)
Gene ID:
8694

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50948-3 targets DGAT1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: DGAT1 Monoclonal antibody, PBS Only (Capture) 60664-5-PBS (1C9D3). 100 µg. Concentration 1 mg/ml.
Detection antibody: DGAT1 Monoclonal antibody, PBS Only (Detector) 60664-6-PBS (2G5H12). 100 µg. Concentration 1 mg/ml.
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50948-3, DGAT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60664-5-PBS. Detection antibody: 60664-6-PBS. Standard:Ag32708. Range: 0.781-100 ng/mL.