

For Research Use Only

OSBPL5 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50959-2

Capture Antibody Information

Catalog Number:

67913-4-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G purification

Clone ID:

2B10E4

Reactivity:

human

Immunogen Catalog Number:

Ag11554

Conjugate:

Unconjugated

Full name:

oxysterol binding protein-like 5

Gene ID:

114879

Detection Antibody Information

Catalog Number:

67913-1-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G purification

Clone ID:

2A1C1

Reactivity:

human, mouse

GenBank:

BC032646

Immunogen Catalog Number:

Ag11554

Conjugate:

Unconjugated

Full name:

oxysterol binding protein-like 5

Gene ID:

114879

Applications

Tested Applications:

Cytometric bead array

Range:

0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50959-2 targets OSBPL5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: OSBPL5 Monoclonal antibody, PBS Only (Capture) 67913-4-PBS (2B10E4). 100 µg. Concentration 1 mg/ml.

Detection antibody: OSBPL5 Monoclonal antibody, PBS Only (Detector) 67913-1-PBS (2A1C1). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50959-2, OSBP.L5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67913-4-PBS. Detection antibody: 67913-1-PBS. Standard: Ag11554. Range: 0.781-100 ng/mL.