

For Research Use Only

FKBP10 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50974-1

Capture Antibody Information

Catalog Number:

68085-2-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G purification

Clone ID:

1F12A11

Reactivity:

human

Immunogen Catalog Number:

Ag31628

Conjugate:

Unconjugated

Full name:

FK506 binding protein 10, 65 kDa

Gene ID:

60681

Detection Antibody Information

Catalog Number:

68085-3-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G Magarose purification

Clone ID:

1E7H11

Reactivity:

human

GenBank:

BC016467

Immunogen Catalog Number:

Ag31628

Conjugate:

Unconjugated

Full name:

FK506 binding protein 10, 65 kDa

Gene ID:

60681

Applications

Tested Applications:

Cytometric bead array

Range:

0.098-25 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50974-1 targets FKBP10 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FKBP10 Monoclonal antibody, PBS Only (Capture) 68085-2-PBS (1F12A11). 100 µg. Concentration 1 mg/mL.

Detection antibody: FKBP10 Monoclonal antibody, PBS Only (Detector) 68085-3-PBS (1E7H11). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50974-1, FKBP10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68085-2-PBS. Detection antibody: 68085-3-PBS. Standard: Ag31628. Range: 0.098-25 ng/mL