

For Research Use Only

PSME1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51008-1

Capture Antibody Information

Catalog Number:
67928-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1E1A5
Reactivity:
human
Immunogen Catalog Number:
Ag30549

Conjugate:
Unconjugated
Full name:
proteasome (prosome, macropain) activator subunit 1 (PA28 alpha)
Gene ID:
5720

Detection Antibody Information

Catalog Number:
67928-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1B4E2
Reactivity:
human
GenBank:
BC007503
Immunogen Catalog Number:
Ag30549

Conjugate:
Unconjugated
Full name:
proteasome (prosome, macropain) activator subunit 1 (PA28 alpha)
Gene ID:
5720

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51008-1 targets PSME1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PSME1 Monoclonal antibody, PBS Only (Capture) 67928-2-PBS (1E1A5). 100 µg. Concentration 1 mg/ml.

Detection antibody: PSME1 Monoclonal antibody, PBS Only (Detector) 67928-3-PBS (1B4E2). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51008-1, PSME1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67928-2-PBS. Detection antibody: 67928-3-PBS. Standard: Ag30549. Range: 0.781-100 ng/mL.