

For Research Use Only

vwf Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51009-6

Capture Antibody Information

Catalog Number:
66682-4-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A purification

Clone ID:
3C12G5

Reactivity:
human

Immunogen Catalog Number:
Ag25578

Conjugate:
Unconjugated

Full name:
von Willebrand factor

Gene ID:
7450

Detection Antibody Information

Catalog Number:
66682-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
4D1C7

Reactivity:
human

GenBank:
Immunogen Catalog Number:
Ag25578

Conjugate:
Unconjugated

Full name:
von Willebrand factor

Gene ID:
7450

Applications

Tested Applications:
Cytometric bead array

Range:
0.195-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51009-6 targets vwf in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: vwf Monoclonal antibody, PBS Only (Capture) 66682-4-PBS (3C12G5). 100 µg. Concentration 1 mg/ml.

Detection antibody: vwf Monoclonal antibody, PBS Only (Detector) 66682-3-PBS (4D1C7). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51009-6, vwf Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66682-4-PBS. Detection antibody: 66682-3-PBS. Standard:Ag25578. Range: 0.195-100 ng/mL.