

For Research Use Only

LTBP2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51013-4

Capture Antibody Information

Catalog Number:
60729-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
1B3C2

Reactivity:
human

Immunogen Catalog Number:
Ag29737

Conjugate:
Unconjugated

Full name:
latent transforming growth factor beta binding protein 2

Gene ID:
4053

Detection Antibody Information

Catalog Number:
60729-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
3D8F3

Reactivity:
human

GenBank:
BC078659

Immunogen Catalog Number:
Ag29737

Conjugate:
Unconjugated

Full name:
latent transforming growth factor beta binding protein 2

Gene ID:
4053

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51013-4 targets LTBP2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LTBP2 Monoclonal antibody, PBS Only (Capture) 60729-1-PBS (1B3C2). 100 µg. Concentration 1 mg/ml.

Detection antibody: LTBP2 Monoclonal antibody, PBS Only (Detector) 60729-2-PBS (3D8F3). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51013-4, LTBP2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60729-1-PBS. Detection antibody: 60729-2-PBS. Standard:Ag29737. Range: 0.391-100 ng/mL.