

For Research Use Only

LGALS9/Galectin-9 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51049-2

Capture Antibody Information

Catalog Number:

60731-1-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G purification

Clone ID:

2G11B11

Reactivity:

human

Conjugate:

Unconjugated

Full name:

lectin, galactoside-binding, soluble, 9

Gene ID:

3965

Detection Antibody Information

Catalog Number:

60731-3-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G Magarose purification

Clone ID:

3D5G6

Reactivity:

human

GenBank:

BC110340

Conjugate:

Unconjugated

Full name:

lectin, galactoside-binding, soluble, 9

Gene ID:

3965

Applications

Tested Applications:

Cytometric bead array

Range:

0.195-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51049-2 targets LGALS9/Galectin-9 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LGALS9/Galectin-9 Monoclonal antibody, PBS Only (Capture) 60731-1-PBS (2G11B11). 100 µg. Concentration 1 mg/ml.

Detection antibody: LGALS9/Galectin-9 Monoclonal antibody, PBS Only (Detector) 60731-3-PBS (3D5G6). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51049-2, LGALS9/Galectin-9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60731-1-PBS. Detection antibody: 60731-3-PBS. Standard:Eg0084. Range: 0.195-100 ng/mL