

For Research Use Only

TIA1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51061-1

Capture Antibody Information

Catalog Number:
68486-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1D9C1

Reactivity:
human

Immunogen Catalog Number:
Ag2778

Conjugate:
Unconjugated

Full name:
TIA1 cytotoxic granule-associated RNA binding protein

Gene ID:
7072

Detection Antibody Information

Catalog Number:
68486-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
2C11F4

Reactivity:
human, mouse, rat

GenBank:
BC015944

Immunogen Catalog Number:
Ag2778

Conjugate:
Unconjugated

Full name:
TIA1 cytotoxic granule-associated RNA binding protein

Gene ID:
7072

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-12.5 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51061-1 targets TIA1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TIA1 Monoclonal antibody, PBS Only (Capture) 68486-2-PBS (1D9C1). 100 µg. Concentration 1 mg/ml.

Detection antibody: TIA1 Monoclonal antibody, PBS Only (Detector) 68486-1-PBS (2C11F4). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51061-1, TIA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68486-2-PBS. Detection antibody: 68486-1-PBS. Standard:Ag2778. Range: 0.098-12.5 ng/mL