

For Research Use Only

# DPEP1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51102-2

## Capture Antibody Information

**Catalog Number:**  
68081-4-PBS

**Host:**  
Mouse

**Isotype:**  
IgG1

**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
1A3B1

**Reactivity:**  
human

**Immunogen Catalog Number:**  
Ag14591

**Conjugate:**  
Unconjugated

**Full name:**  
dipeptidase 1 (renal)

**Gene ID:**  
1800

## Detection Antibody Information

**Catalog Number:**  
68081-5-PBS

**Host:**  
Mouse

**Isotype:**  
IgG1

**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
1F7C10

**Reactivity:**  
human

**GenBank:**  
BC017023

**Immunogen Catalog Number:**  
Ag14591

**Conjugate:**  
Unconjugated

**Full name:**  
dipeptidase 1 (renal)

**Gene ID:**  
1800

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
1.563-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51102-2 targets DPEP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DPEP1 Monoclonal antibody, PBS Only (Capture) 68081-4-PBS (1A3B1). 100 µg. Concentration 1 mg/ml.

Detection antibody: DPEP1 Monoclonal antibody, PBS Only (Detector) 68081-5-PBS (1F7C10). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51102-2, DPEP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68081-4-PBS. Detection antibody: 68081-5-PBS. Standard:Ag14591. Range: 1.563-100 ng/mL.