

For Research Use Only

p130Cas / BCAR1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51156-1

Capture Antibody Information

Catalog Number:
67215-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1E7D8
Reactivity:
human
Immunogen Catalog Number:
Ag10604

Conjugate:
Unconjugated
Full name:
breast cancer anti-estrogen resistance 1
Gene ID:
9564

Detection Antibody Information

Catalog Number:
67215-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1B9G5
Reactivity:
human
GenBank:
BC062556
Immunogen Catalog Number:
Ag10604

Conjugate:
Unconjugated
Full name:
breast cancer anti-estrogen resistance 1
Gene ID:
9564

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-12.5 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51156-1 targets p130Cas / BCAR1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: p130Cas / BCAR1 Monoclonal antibody, PBS Only (Capture) 67215-2-PBS (1E7D8). 100 µg. Concentration 1 mg/mL.

Detection antibody: p130Cas / BCAR1 Monoclonal antibody, PBS Only (Detector) 67215-3-PBS (1B9G5). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51156-1, p130Cas / BCAR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67215-2-PBS. Detection antibody: 67215-3-PBS. Standard:Ag10604. Range: 0.781-12.5 ng/mL