

For Research Use Only

C5AR1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51210-4

Capture Antibody Information

Catalog Number:
60761-6-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
1F2G1

Reactivity:
human

Immunogen Catalog Number:
Ag29140

Conjugate:
Unconjugated

Full name:
complement component 5a receptor 1

Gene ID:
728

Detection Antibody Information

Catalog Number:
60761-2-PBS

Host:
Mouse

Isotype:
IgG2b

Purification Method:
Protein A Magarose purification

Clone ID:
1G5F10

Reactivity:
human

GenBank:
BC008982

Immunogen Catalog Number:
Ag29140

Conjugate:
Unconjugated

Full name:
complement component 5a receptor 1

Gene ID:
728

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51210-4 targets C5AR1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: C5AR1 Monoclonal antibody, PBS Only (Capture) 60761-6-PBS (1F2G1). 100 µg. Concentration 1 mg/mL.

Detection antibody: C5AR1 Monoclonal antibody, PBS Only (Detector) 60761-2-PBS (1G5F10). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51210-4, CSAR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60761-6-PBS. Detection antibody: 60761-2-PBS. Standard:Ag29140. Range: 0.781-100 ng/mL.