

For Research Use Only

# FREQ Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51282-1

## Capture Antibody Information

Catalog Number:  
67616-1-PBS

Host:  
Mouse

Isotype:  
IgG2a

Purification Method:  
Protein A purification

Clone ID:  
2F10C8

Reactivity:  
human, mouse, rat, pig

Immunogen Catalog Number:  
Ag29808

Conjugate:  
Unconjugated

Full name:  
frequenin homolog (Drosophila)

Gene ID:  
23413

## Detection Antibody Information

Catalog Number:  
67616-2-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
2D7G6

Reactivity:  
human

GenBank:  
BC004856

Immunogen Catalog Number:  
Ag29808

Conjugate:  
Unconjugated

Full name:  
frequenin homolog (Drosophila)

Gene ID:  
23413

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51282-1 targets FREQ in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NCS1 Monoclonal antibody, PBS Only (Capture) 67616-1-PBS (2F10C8). 100 µg. Concentration 1 mg/mL.

Detection antibody: FREQ Monoclonal antibody, PBS Only (Detector) 67616-2-PBS (2D7G6). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP51282-1, FREQ Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67616-1-PBS. Detection antibody: 67616-2-PBS. Standard:Ag29808. Range: 1.563-100 ng/mL.