

IL-1 beta Hybrid Matched Antibody Pair, PBS Only

Catalog Number:MP80001-6

Capture Antibody Information

Catalog Number:

82696-8-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

2G9

Reactivity:

Human

Conjugate:

Unconjugated

Full name:

interleukin 1, beta

Gene ID:

3553

Detection Antibody Information

Catalog Number:

66737-1-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein A purification

Clone ID:

2A1B4

Reactivity:

Human

GenBank:

BC008678

Immunogen Catalog Number:

Ag26030

Conjugate:

Unconjugated

Full name:

interleukin 1, beta

Gene ID:

3553

Applications

Tested Applications:

Sandwich ELISA

Range:

0.098-12.5 ng/mL (Sandwich ELISA)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP80001-6 targets IL-1 beta in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: IL1B Recombinant antibody, PBS Only (Capture) 82696-8-PBS (2G9). 100 µg. Concentration 1 mg/ml.

Detection antibody: IL-1 Beta Monoclonal antibody, PBS Only (Detector) 66737-1-PBS (2A1B4). 100 µg. Concentration 1 mg/ml.

Alternative IL-1 beta matched antibody pairs: MP00331-1, MP50134-1, MP80001-1, MP80001-2, MP80001-3, MP80001-4, MP80001-5

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP80001-6, IL-1 beta Hybrid Matched Antibody Pair, PBS Only. Capture antibody: 82696-8-RR. Detection antibody: HRP-conjugated 66737-1-PBS. Standard:HZ-1164. Range: 0.098-12.5 ng/mL