

# BACH1 Hybrid Matched Antibody Pair, PBS Only

Catalog Number:MP80007-3

## Capture Antibody Information

**Catalog Number:**

66762-1-PBS

**Host:**

Mouse

**Isotype:**

IgG1

**Purification Method:**

Protein A purification

**Clone ID:**

2D12A4

**Reactivity:**

Human, Rat, mouse, pig, rabbit

**Immunogen Catalog Number:**

Ag6011

**Conjugate:**

Unconjugated

**Full name:**

BTB and CNC homology 1, basic leucine zipper transcription factor 1

**Gene ID:**

571

## Detection Antibody Information

**Catalog Number:**

82842-2-PBS

**Host:**

Rabbit

**Isotype:**

IgG

**Purification Method:**

Protein A purification

**Clone ID:**

4O2

**Reactivity:**

Human

**GenBank:**

BC063307

**Immunogen Catalog Number:**

Ag5128

**Conjugate:**

Unconjugated

**Full name:**

BTB and CNC homology 1, basic leucine zipper transcription factor 1

**Gene ID:**

571

## Applications

**Tested Applications:**

Cytometric bead array

**Range:**

0.313-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP80007-3 targets BACH1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BACH1 Monoclonal antibody, PBS Only (Capture) 66762-1-PBS (2D12A4). 100 µg. Concentration 1 mg/ml.

Detection antibody: BACH1 Recombinant antibody, PBS Only (Detector) 82842-2-PBS (4O2). 100 µg. Concentration 1 mg/ml.

Alternative BACH1 matched antibody pairs: MP80007-1, MP80007-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology.

Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**

Store at -80°C.

**Storage buffer:**

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP80007-3, BACH1 Hybrid Matched Antibody Pair, PBS Only. Capture antibody: 66762-1-PBS. Detection antibody: 82842-2-PBS. Standard: Ag6011. Range: 0.313-40 ng/mL