

For Research Use Only

PE Anti-Mouse CD48 (HM48-1)

Catalog Number: PE-65227



Basic Information

Catalog Number:

PE-65227

Size:

100ug, 0.2 mg/ml

Source:

Armenian Hamster

Isotype:

IgG

GenBank Accession Number:

BC060977

GeneID (NCBI):

12506

UNIPROT ID:

P18181

Full Name:

CD48 antigen

Purification Method:

Affinity purification

CloneNo.:

HM48-1

Excitation/Emission maxima
wavelengths:

496 nm, 565 nm / 578 nm

Applications

Tested Applications:

FC

Species Specificity:

mouse

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ C57BL/6 mouse splenocytes were surface stained with 0.1 ug PE Anti-Mouse CD48 (PE-65227, Clone: HM48-1), or 0.1 ug PE Armenian Hamster IgG Isotype Control (PIP) (PE-65210, Clone: PIP). Cells were not fixed.