

For Research Use Only

PE Anti-Rat CD43 (W3/13) Mouse IgG2a Recombinant Antibody

Catalog Number: PE-65607



Basic Information

Catalog Number:

PE-65607

Size:

100ug , 200 µg/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

GeneID (NCBI):

24796

Full Name:

sialophorin

Calculated MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

W3/13

Excitation/Emission maxima wavelengths:

496 nm, 565 nm / 578 nm

Applications

Tested Applications:

FC

Species Specificity:

rat

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ rat splenocytes were surface stained with 0.1 ug PE Anti-Rat CD43 (W3/13) Mouse IgG2a Recombinant Antibody (PE-65607, Clone: W3/13) or 0.1 ug PE Mouse IgG2a Isotype Control (C1.18.4) (PE-65208, Clone: C1.18.4). Cells were not fixed.



1x10⁶ rat splenocytes were surface stained with 0.1 ug PE Anti-Rat CD43 (W3/13) Mouse IgG2a Recombinant Antibody (PE-65607, Clone: W3/13) (red), or 0.1 ug PE Mouse IgG2a Isotype Control (C1.18.4) (PE-65208, Clone: C1.18.4) (blue). Cells were not fixed.