

For Research Use Only

FcZero-rAb™ PE Anti-Human CD70 Rabbit Recombinant Antibody

Catalog Number: PE-FcA98098



Basic Information

Catalog Number:

PE-FcA98098

Size:

100tests , 5 ul/test

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC000725

GeneID (NCBI):

970

Full Name:

CD70 molecule

Calculated MW:

21 kDa

Purification Method:

Protein A purification

CloneNo.:

241325D8

Excitation/Emission maxima
wavelengths:

496 nm, 565 nm / 578 nm

Applications

Tested Applications:

FC

Species Specificity:

human

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ U266 cells were surface stained with 5 ul PE Anti-Human CD70 Rabbit Recombinant Antibody (PE-FcA98098, Clone: 241325D8) (red), or 5 ul PE Rabbit IgG Isotype Control Recombinant Antibody (PE-FcA98136, Clone: 240953C9) (blue). Cells were not fixed.