

For Research Use Only

# FcZero-rAb™ PE Anti-Human CD89 Rabbit Recombinant Antibody

Catalog Number: PE-FcA98105



## Basic Information

<b>Catalog Number:</b> PE-FcA98105	<b>GenBank Accession Number:</b> BC027953	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100tests, 5 ul/test	<b>GeneID (NCBI):</b> 2204	<b>CloneNo.:</b> 241426E3
<b>Source:</b> Rabbit	<b>Full Name:</b> Fc fragment of IgA, receptor for	<b>Excitation/Emission maxima wavelengths:</b> 496 nm, 565 nm / 578 nm
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 287 aa, 32 kDa	

## Applications

**Tested Applications:**  
FC

**Species Specificity:**  
human

## Storage


**Storage:**  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)     E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)     W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> human peripheral blood leukocytes were surface stained with 5 ul PE Anti-Human CD89 Rabbit Recombinant Antibody (PE-FcA98105, Clone: 241426E3), or 5 ul PE Rabbit IgG Isotype Control Recombinant Antibody (PE-FcA98136, Clone: 240953C9). Cells were not fixed.