

For Research Use Only

Foxp3 / Transcription Factor Fix/Perm Diluent (1X)

Catalog Number: PF00011-B



Description	Foxp3 / Transcription Factor Fix/Perm Diluent (1X) Intended for use as a diluent for Foxp3 / Transcription Factor Fix/Perm Concentrate (4X)
Package	160 mL
Storage	4°C
Expiration Date	6 months after receiving.

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB(1-888-478-4522)(toll free in USA), or
1(312) 455-8498(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.