

For Research Use Only

RBC Lysis Buffer (10X)

Catalog Number: PF00013



www.ptglab.com

Description

RBC Lysis Buffer (10X) is a concentrated ammonium chloride-based lysing reagent. The diluted 1X working solution will lyse red blood cells in single cell suspensions with minimal effects on leukocytes. RBC Lysis Buffer (10X) does not contain a fixative so the cells remain viable after.

Package

100 mL

Storage

4°C

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB(1-888-478-4522)(toll free in USA), or
1(312) 455-8498(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.