

For Research Use Only

Spot-Label® ATTO488



www.ptglab.com

Catalog Number: eba488

1 Publications

Basic Information

Catalog Number:
eba488

Applications:
IF

Host:
Alpaca

Conjugate:
ATTO 488

Type :

Nanobody

Class:

Recombinant

RRID:

AB_2827569

Description

anti-Spot-Tag® VHH/ Nanobody conjugated to fluorescent dye for immunofluorescence, microscopy, and immunoblotting of Spot-tagged proteins

Specificity/Target

Spot-Tag® (PDRVRAVSHWSS)

Physical State

Liquid

Suggested Dilution

IF: 1:1,000

Affinity (K_D)

Dissociation constant K_D of 6 nM

Storage

Storage:

Shipped at ambient temperature. Upon receipt store at +4°C. Stable for 6 month. Do not freeze. Protect from light.

Storage Buffer:

1x PBS, 0.09% sodium azide

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: proteintech@ptglab.com
1(312) 455-8498 (outside USA) W: www.ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

	Spot-Trap	Spot-Cap	Spot-Label	Spot V _H H
Immunoprecipitation (IP) / Co-IP / Co-IP/MS	✓			
Protein purification		✓		
Western Blotting			✓	✓
Immunofluorescence			✓	✓

SPOT tag related products