

For Research Use Only

# Spot-Label® Alexa Fluor® 568



www.ptglab.com

Catalog Number: **ebAF568**

4 Publications

## Basic Information

**Catalog Number:**  
ebAF568

**Applications:**  
IF

**Host:**  
Alpaca

**Conjugate:**  
Alexa Fluor® 568

**Type :**

Nanobody

**Class:**

Recombinant

**RRID:**

AB\_2889375

## Description

anti-Spot-Tag VHH/ Nanobody conjugated to fluorescent dye for immunofluorescence, microscopy, and immunoblotting of Spot-tagged proteins

## Specificity/Target

Spot-Tag® (PDRVRAVSHWSS)

## Physical State

Liquid

## Suggested Dilution

IF: 1:800

## DOL

1 fluorophore per bivalent Nanobody, site-directed conjugation

## Affinity (K<sub>D</sub>)

Dissociation constant K<sub>D</sub> of 6 nM

## Storage

**Storage:**

Shipped at ambient temperature. Store at -20°C/-4°F. Stable for 6 month. Protect from light.

**Storage Buffer:**

PBS, 5 mM EDTA, preservative: 0.09 % sodium azide

For technical support and original validation data for this product please contact


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
1(312) 455-8498 (outside USA) W: [www.ptglab.com](http://www.ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

	Spot-Trap	Spot-Cap	Spot-Label	Spot V <sub>H</sub> H
Immunoprecipitation (IP) / Co-IP / Co-IP/MS	✓			
Protein purification		✓		
Western Blotting			✓	✓
Immunofluorescence			✓	✓

SPOT tag related products



Left: HeLa cells transiently expressing pSpot Spot-Tag<sup>®</sup>-Actin (ev-31) were immunostained with Spot-Label<sup>®</sup> Alexa Fluor<sup>®</sup> 568 (red, ebAF568, 1:800) and DAPI (blue). Right: Control staining with Spot-Label<sup>®</sup> Alexa Fluor<sup>®</sup> 568 (red, ebAF568, 1:800) and DAPI (blue) of untransfected HeLa cells showing no unspecific background. Scale bar, 10  $\mu$ m. Images were acquired with the Thermo Scientific CellInsight CX7, 20X objective.