

Catalog Number: **hta**

Basic Information

Catalog Number:

hta

Applications:

IP, CoIP

Conjugate:

Agarose beads; bead size: ~ 90 µm (cross-linked 4 % agarose beads)

Host:

Alpaca

Type:

Nanobody

Class:

Recombinant

Description

Mdm4/HdmX-Trap Agarose is an affinity resin for IP of Mdm4/HdmX-fusion proteins.

Specificity/Target

Double minute 4 protein (Mdm4/HdmX), human

Note that binding of Mdm4/HdmX-Trap interferes with binding of p53 to Mdm4 HdmX.

Elution buffer

SDS sample buffer

0.2 M glycine pH 2.5

Storage

Storage:

Upon receipt store at +4°C. Do not freeze!

Storage Buffer:

20% ethanol

Selected Validation Data

