

For Research Use Only

Mdm4/ HdmX-Trap Agarose, Kit for Immunoprecipitation



www.ptglab.com

Catalog Number: **htak**

Basic Information

Catalog Number:

htak

Applications:

IP, CoIP

Conjugate:

Agarose beads; bead size: ~ 90 µm (cross-linked 4 % agarose beads)

Host:

Alpaca

Type:

Nanobody

Class:

Recombinant

Description

Mdm4/HdmX-Trap Agarose is an affinity resin for IP of Mdm4/HdmX-fusion proteins.

Specificity/Target

Double minute 4 protein (Mdm4/HdmX), human

Note that binding of Mdm4/HdmX-Trap interferes with binding of p53 to Mdm4 HdmX.

Elution buffer

SDS sample buffer

0.2 M glycine pH 2.5

Storage

Storage:

Upon receipt store at +4°C. Do not freeze!

Storage Buffer:

20% ethanol

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: proteintech@ptglab.com
1(312) 455-8498 (outside USA) W: www.ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

