

For Research Use Only

mNeonGreen RecAb, VHH-rabbit IgG Fc fusion



www.ptglab.com

Catalog Number: nfrb

1 Publications

Basic Information

Catalog Number:
nfrb

Clone No.:
CTK0203

Applications:
IF

Host/IsoType:
Rabbit/IgG

Conjugate:
Unconjugated

Type:

Primary antibody; chimeric fusion of Nanobody (VHH) and rabbit IgG Fc domain

Class:

Recombinant

RRID:

AB_2889369

Purification Method:

Recombinantly expressed and purified by IMAC

Description

Rabbit monoclonal antibody to mNeonGreen. Recombinant fusion of anti-mNeonGreen VHH [CTK0203] and rabbit IgG Fc domain.

Specificity/Target

mNeonGreen fluorescent protein

Recommended Dilution

IF/ICC: 1:200 - 1:2,000

Storage

Storage:

Shipped on ice. Aliquot upon receipt and store at -20°C/-4°F. Avoid freeze-thaw cycles.

Storage Buffer:


PBS, 10% glycerol, Preservative: 0.02 % sodium azide, safety datasheet (SDS): sodium azide

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: proteintech@ptglab.com
1(312) 455-8498 (outside USA) W: www.ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells transiently expressing Lamin-Chromobody-mNeonGreen were immunostained with mNeonGreen recombinant antibody, VHH-rabbit IgG Fc fusion [CTK0203] (nfrb, 1:1,000) and Nano-Secondary alpaca anti-human IgG/anti-rabbit IgG, recombinant VHH, Alexa Fluor 647 [CTK0101, CTK0102] (srbAF647-1, 1:1,000). Scale bar, 10 μ M.