

For Research Use Only

p53 N-term-Trap Agarose, Kit for Immunoprecipitation



www.ptglab.com

Catalog Number: **ptak**

Basic Information

Catalog Number:

ptak

Applications:

IP, CoIP

Conjugate:

Agarose beads; bead size: ~ 90 µm (cross-linked 4 % agarose beads)

Host:

Alpaca

Type:

Nanobody

Class:

Recombinant

Description

p53-N-term-Trap Agarose is an affinity reagent for IP of p53. It consists of a p53 Nanobody/ VHH coupled to agarose beads, and which binds to the N-terminus of p53.

Specificity/Target

p53 (isoforms: alpha, beta, and gamma); epitope 1-81 aa

Elution buffer

SDS sample buffer

0.2 M glycine pH 2.5

Storage

Storage:

Upon receipt store at +4°C. Do not freeze!

Storage Buffer:


20% ethanol

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: proteintech@ptglab.com
1(312) 455-8498 (outside USA) W: www.ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



p53-N-term-Trap for immunoprecipitation of p53 isoform alpha, beta and gamma. I: Input, FT: Flow-through, B: Bound.