

chromotek® Strep-Peptide

Product Code: qp

Product Information

Description: HPLC-purified Strep-peptide

Sequence: SAWSHPQFEK

Molecular Weight: 1216 Da

Form: HPLC-purified peptide, lyophilized powder

Purity: 95 %

Application: Elution of proteins fused to SAWSHPQFEK - (also known as Strep-Tag®) bound to Strep-NanoTrap.

Reconstitution: Reconstitute peptide in 8.2 ml PBS, which results in a final concentration of 100 µM. Vortex for 1 min to dissolve the powder.

Shipment: Shipped at ambient temperature

Storage Condition: Upon receipt store at -20°C. Aliquot upon reconstitution and store at -20°C. Avoid repeated freeze/thaw cycles.

Stability: Stable for 1 year upon receipt

Product Sizes

Product	Product Code	Size
Strep-Peptide	qp-1	1 mg

Related Products

Product	Code
Strep-NanoTrap Agarose	qta
Strep-NanoTrap Agarose Kit	qtak

chromotek® Strep-Peptide

Product Code: qp

Contact

Proteintech North America (HQ)

Proteintech Group, Inc.
5500 Pearl Street, Suite 400
Rosemont, IL 60018 USA

1-888-472-4522
proteintech@ptglab.com

Proteintech Europe

Transmission (6th Fl)
6 Atherton Street
M3 3GS, Manchester, UK

+44 161 839 3007
europe@ptglab.com

ChromoTek & Proteintech Germany

Fraunhoferstr. 1
82152, Planegg-Martinsreid
Germany

+49 89 124 148 850
germany@ptglab.com

Disclaimer

For research use only (RUO), not for diagnostic or therapeutic use.

ChromoTek® is a registered trademark of ChromoTek GmbH, part of Proteintech Group.
Strep-tag® and Twin-Strep-Tag® are registered trademarks of IBA Lifesciences GmbH.