

For Research Use Only

Nano-Secondary® anti-mouse IgG1, recombinant VHH, for 2x Cys conjugation [CTK0104]

Catalog Number: smsG1Cys2-2



www.ptglab.com

Basic Information

Catalog Number:

smsG1Cys2-2

Applications:

IF, WB, FC

Host:

Alpaca

Conjugate:

N- and C-terminal cysteine conjugation with thiol-reactive reagents

Type:

Nano-Secondary®, secondary Nanobody (VHH)

RRID:

AB_2864264

Description

Monovalent, recombinant secondary single domain antibody to mouse IgG1: alpaca monoclonal Nanobody, Fc-specific, for cysteine conjugation with thiol-reactive reagents, e.g. maleimides

Specificity/Target

Fc-fragment of mouse IgG1

Affinity

K_D = 0.6 nM

Storage

Storage:

Shipped on dry ice. Store at -20°C/-4°F. Avoid freeze-thaw cycles.

Storage Buffer:

10 mM HEPES pH 7.0, 500 mM NaCl, 1 mM TCEP
Preservative: 0.09 % Sodium azide

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: proteintech@ptglab.com
1(312) 455-8498 (outside USA) W: www.ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.