

For Research Use Only

# anti-V5 recombinant VHH, unconjugated



[www.ptglab.com](http://www.ptglab.com)

Catalog Number: v5t

## Basic Information

**Catalog Number:**  
v5t

**Applications:**  
Suitable for conjugation of dyes, biotin, beads, surfaces etc. via NHS ester reaction

**Host:**  
Alpaca

**Conjugate:**  
Unconjugated

**Type :**  
Nanobody

**Class:**  
Recombinant

**RRID:**  
AB\_2868499

**Molecular Weight:**  
14.6 kDa

## Description

Alpaca anti-V5 VHH, purified recombinant binding protein

## Specificity/Target

V5-tag sequence GKIPNPPLGLDST

## Affinity ( $K_D$ )

Dissociation constant  $K_D$  of 40 nM

## Storage

**Storage:**  
Shipped at ambient temperature . Upon receipt store at 4°C; stable for one year.

**Storage Buffer:**  
25mM TAPS pH 8.5, 137 mM NaCl, 5 mM EDTA, preservative: 0.09 % sodium azide

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
1(312) 455-8498 (outside USA) W: [www.ptglab.com](http://www.ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Structure of VHH.